



## Conformity Statement

(1)

(2) Equipment and protective systems intended for use in potentially explosive atmospheres - **Directive 2014/34/EU**

(3) Certificate number: **SEV 08 ATEX 0150 X**

(4) Product: Freezer LIEBHERR  
Type 4 6183 4 (LGUex 1500)

(5) Manufacturer: Liebherr-Hausgeräte Lienz GmbH

(6) Address: Dr. Hans-Liebherr-Strasse 1, 9900 Lienz, AUSTRIA

(7) This product and any acceptable variation thereto are specified in the schedule to this certificate and the documents therein referred to.

(8) Electrosuisse SEV certifies that this product has been found to comply with the essential health and safety requirements relating to the design and construction of equipment and protective systems intended for use in potentially explosive atmospheres given in Annex II to Directive 2014/34/EU of the European parliament and of the Council, dated 26 February 2014.

The examination and test results are recorded in confidential report no. 08-IK-0329.01 + Extension1

(9) Compliance with the essential health and safety requirements has been assured by compliance with:

**EN 60079-0:12 + A11:13      EN 60079-15:10**

Except in respect of those requirements listed at item 18 of the schedule.

(10) If the sign «X» is placed after the certificate number, it indicates that the product is subjected to special conditions for safe use specified in the schedule to this certificate.

(11) This Conformity Statement relates only to the design and construction of the specified product and not to specific items of product subsequently manufactured.

(12) The marking of the product shall include the following:



**II 3G, Ex nA II T6 Gc**

Electrosuisse  
Notified Body ATEX

Martin Plüss  
Product Certification



(13) **Appendix**

(14) **Conformity Statement**

(15) **Description of product**

The freezer LIEBHERR Type 4 6183 4 ( LGUex 1500 ) is to installed in non hazardous areas. The interior room is free of ignition sources and therefore intended for storage of combustible materials, partially in open containers, and is classified as explosive atmosphere of zone 2. All electrical components, apart from a air sensor, are not installed in the hazardous areas.

Rating(s):

AC 220 -240 V, 50 Hz, max. 1 A  
AC 100 V, 50-60 Hz, max. 2 A

Remark:

With start-up of the combined refrigerator/freezer the notes in the instructions of the manufacturer are to be observed.

(16) **Report number** 08-IK-0329.01 + Extension1

(17) **Specific conditions of use**

WARNING – Potential electrostatic charging hazard – See instruction manual!

(18) **Essential health and safety requirements**

In addition to the essential health and safety requirements (EHSRs) covered by the standards listed at item 9, the following are considered relevant to this product, and conformity is demonstrated in the report:

Clause	Subject
None	

(19) **Drawings and documents**

Number	Sheet	Issue date	Description
See test report "Manufacturer's Documents"			